

AutoEMage user manual

CHENG Yuanhao

Version 1.0.0

Contact: chengyuanhao@iphy.ac.cn

September 6, 2023

Contents

1	Introduction	2
2	System Requirements and Installation	2
3	Usage	3
3.1	Graphical User Interface and Functions	3
3.1.1	Registration and Login	3
3.1.2	File Transfer	3
3.1.3	Variables Display	3
3.1.4	Autopicking and 2D class averaging	4
3.1.5	Class ranking and 3D model generation	5
3.1.6	Outliers Detection and Email Reminder	5
3.2	Incomplete Functions	5

1 Introduction

AutoEMage is a software that automates file transfer, motion correction, CTF estimation, image display, outlier detection, particle picking, 2D class averaging, class ranking and 3D model generation during cryo-EM data acquisition. Through specific settings, file transfer, data preprocessing, real-time display and monitoring can be executed simultaneously with data acquisition. In this process, AutoEMage does not require manual intervention or surveillance, which lives up to its name. Once outliers are detected, it will send emails to users so that they can come back in time and have a check. Moreover, it provides users with a graphical user interface, in order that they can not only see real-time preprocessed micrographs, but also some results of data processing.

2 System Requirements and Installation

AutoEMage is developed using Python3 and Perl on Linux, including a few programs such as MotionCor2, CTFFIND4, UCSF ChimeraX and some programs of RELION and IMOD, while its graphical user interface is based on PyQt6. Therefore, these programming languages and cryo-EM programs should be installed before using AutoEMage, which can be downloaded through the following links:

- PyQt6: <https://pypi.org/project/PyQt6/> Note that PyQt6 may be successfully installed by inputting “python3 -m pip install PyQt6” on the command line.
- Python3: <https://www.python.org/>
- Perl: <https://www.perl.org/>
- MotionCor2: <https://hpc.nih.gov/apps/MotionCor2.html>
- CTFFIND4: <https://grigoriefflab.umassmed.edu/ctffind4>
- RELION: <https://relion.readthedocs.io/en/latest/Installation.html>
- IMOD: <https://bio3d.colorado.edu/imod/download.html>
- UCSF ChimeraX: <https://www.cgl.ucsf.edu/chimerax/>

Warning: MotionCor2 supports cuda of different versions, but RELION does not support cuda 12.1.

In addition, AutoEMage can be installed in a virtual environment created by “virtualenv”. Some libraries such as PyQt6-WebEngine, matplotlib, scipy, tiffle, scikit-image and mrcfile need manual installation. To install the aforementioned libraries, users can enter the following command in the terminal: python3 -m pip install PyQt6-WebEngine matplotlib mrcfile tiffle scipy scikit-image. After successful installation, directories of MotionCor2, CTFFIND4 and AutoEMage can be added to environment variables in this way: add the following lines to the .bashrc file:

```
#add MotionCor2 path
PATH=$PATH:/work/Softwares/MotionCor2_1.6.4_Mar31_2023
#add CTFFIND 4 path
PATH=$PATH:/work/Softwares/ctffind-4.1.14-linux64/bin
#add AutoEMage path
PATH=$PATH:/work/Softwares/autoemage_2023
# add IMOD
export IMOD_DIR=/usr/local/imod_4.11.24
if [ -e $IMOD_DIR/IMOD-linux.sh ] ; then source $IMOD_DIR/IMOD-linux.sh ; fi
# Setup |RELION| if not already done so
if [ "" == "echo $PATH | grep /work/Documents/relion/build/bin" ]; then
PATH=$PATH:/work/Documents/relion/build/bin
fi
if [ "" != "echo $LD_LIBRARY_PATH" ]; then
if [ "" == "echo $LD_LIBRARY_PATH | grep /work/Documents/relion/build/lib" ]; then
LD_LIBRARY_PATH=$LD_LIBRARY_PATH:/work/Documents/relion/build/lib
fi
else
```

```
export LD_LIBRARY_PATH=$LD_LIBRARY_PATH:/work/Documents/relion/build/lib
fi
```

3 Usage

3.1 Graphical User Interface and Functions

3.1.1 Registration and Login

Users can run “autoemage.py” on the command line to open the user graphical interface. Up to now, the login interface is only used for keeping a record of users’ emails, in order that AutoEMage can send emails and store data for each user respectively. Sending email alerts is one convenient function of AutoEMage whose email address is cryo_EM_2022@163.com. To make use of this function, users had better type in their current email addresses during registration. The login interface is shown in Figure 1.



Figure 1: AutoEMage login interface

Following successful login, the main interface shows up (Figure 2). During data acquisition, the left panel will display several physical quantities varying with time and the right panel will exhibit plots of “preprocessed micrograph–CTF fit–global drift path–patch frame drift”.

3.1.2 File Transfer

In the main interface, users can click the “File Transfer” button at the top left corner to open the window to input some necessary experimental values, as shown in Figure 3. The EM directory refers to the directory where the electron microscope automatically saves images, and the destination directory is the directory of users’ mobile disk drive or any other directory to store data. The job name should contain no space. After inputting the parameters for the first time, users can save the current profile to a text file for usage in the future. After inputting necessary values, users can start automatic file transfer by clicking “OK”, while the progress will update in the progress bar. As the file transfer goes on, gain correction, motion correction and CTF calculation will simultaneously begin and the results will show up once finished, as shown in Figure 4.

3.1.3 Variables Display

The left panel displays a few physical variables changing with time during data acquisition. Time sliders below the graph are used for displaying variables in a specific time slot. As thousands of images are usually collected for a project, thousands of data points will be produced. If all of them are plotted on a single graph, the graph will be too crowded to exhibit detailed variations. The first time slider controls the start of the time slot, and the second slider determines the end of the time slot.




Figure 2: The main GUI of AutoEMage




Figure 3: The input UI of experimental parameters

3.1.4 Autopicking and 2D class averaging

To utilize the autopicking and 2D class averaging functions, users simply need to input the smallest and largest diameter of the target particles in the provided interface. These functions rely on programs from RELION-4.0, such as `relion_autopick`, `relion_preprocess`, and `relion_refine`. Autopicking is executed twice for each micrograph. Initially, laplacian of Gaussian picking is employed, and the total number of particles will be displayed at the bottom right corner of the main GUI. If the cumulative number of particles from multiple




Figure 4: The main GUI of AutoEMage during data acquisition

micrographs exceeds a predefined value (e.g., 20,000, adjustable), 2D classification is performed and the results are used as references for subsequent autopicking. For every 20,000 particles, reference-based autopicking is carried out, followed by 2D classification. Users familiar with RELION can customize the corresponding script ([autoimage_threads.py](#)) according to their specific samples to optimize the performance of these programs.

3.1.5 Class ranking and 3D model generation

After the completion of 2D classification, an automatic class ranking is carried out using `reliion_class_ranker`. This ranking assigns a score to each class, and the "good" classes are selected to generate a 3D initial model using `reliion_refine`. For particles beyond the initial 20,000, 3D model generation occurs every 20,000 particles, with the results visualized in UCSF ChimeraX. Users have the flexibility to adjust the default processing parameters in the script [autoimage_threads.py](#).

3.1.6 Outliers Detection and Email Reminder

Data collection usually stores thousands of images, so outliers are inevitable. When outliers appear frequently within a short amount of time, for instance 10 out of 50 images, AutoEMage will send an email to users so that they can come back to check the EM settings as well as their samples, as shown in Figure 5.

Up to now, AutoEMage can detect two kinds of outliers: images with large drift and images with low CTF scores, as shown in Figure 6. These images are unlikely to contribute to the high-resolution reconstruction.

When the task is finished or the hard disk drive is nearly full, AutoEMage will also send a reminding email to users.

3.2 Incomplete Functions

AutoEMage is an open source software for academic users, so its biggest advantage is the freedom that users have to modify and add any function they need. As this software was used by different groups in their own labs, they have left their special scripts that have witnessed the development of this software as well as the ever growing requirements of a practical software in the field of cryo-EM. Here I point out these special codes for users' reference. For example, lines 242-340 in the script [Auto_ctf_find.pl](#) and the script [Stigma_check_hxj_v1.pl](#). These are also reflected in "EM settings" located in the toolbar of the main interface.

Outliers warning

cryo_EM_2022@163.com 发送给 程渊溟 包含 10 个附件

Dear scientists,

Low-score data shows up after CtfFind. Please examine EM settings or your samples in time.

AutoEMage

普通附件(10 个) 全部下载 全部保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0014_all.png (187.17K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0026_all.png (201.31K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0029_all.png (191.64K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0031_all.png (195.54K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0034_all.png (228.53K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0036_all.png (195.61K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0042_all.png (222.33K)

下载

预览

保存到个人网盘



/work2/cyh/Pictures/Data/Outliers/CtfFind/Data_0047_all.png (204.18K)

Figure 5: Outliers email




Figure 6: Two kinds of outliers. On the left is an image with large drift. On the right is an image with a low CTF score.